

Modified McMaster Egg Count Procedure

Principle

Many eggs and larvae will float in certain sugar/salt flotation solutions because they have a lower specific gravity than the solutions they are floating in. This procedure employs a special counting chamber that allows you to count the number of eggs present in .300 ml of available solution in the chambers.

Materials and Equipment

- 1) Tongue depressors
- 2) Plastic cups
- 3) Cheesecloth (5 X 5 double layer)
- 4) Sugar/Salt flotation solution (1125 g sugar, 748 g salt, 2 L distilled water)
- 5) Pipettes
- 6) McMaster egg counting chamber slides
- 7) Microscope (10X objective)
- 8) Scale
- 9) Counter

Procedure

- 1) Weigh out 4 g fecal material (approximately 1 tbsp).
- 2) To one plastic cup, add weighed sample and 26 ml flotation solution.
- 3) Thoroughly break up fecal sample, mix, and strain through cheesecloth.
- 4) Immediately fill both sides of the McMasters chamber.
 - a. Fill pipette from the surface since the theory is the eggs will float to the top.
 - b. It may also help to get a better sample size by emptying anything that is left in the pipette after the 1st chamber is filled, waiting a few seconds (approximately 30), and then refilling the pipette again from the surface to fill the second side of the chamber.
- 5) Let sample sit in the chamber for 2-5 minutes to allow eggs to float to the surface of the chamber.
- 6) Using the 10X objective (lowest objective on microscope) count the number of eggs in the sample going up and down the 6 lanes inscribed on each side of the counting chamber.
 - a. It will help to focus on and begin counting at the corner square of your choice.
- 7) Multiply the number of eggs counted by 25.

Calculation Explanation

- 1) 4 g fecal material + 26 ml sugar/salt solution = 30 ml
- 2) Amount of sample examined in McMaster slide = .300 ml (0.150 ml/slide)
- 3)

30 (total volume)		100
-----	=	-----
.300 (chamber)		4 (g fecal material)
	=	25 (multiplier)

Source: Dr. Dahlia O'Brien, Delaware State University